

Tailor™ Flexible Annuloplasty Ring

Flexible Rings



Product Highlights

- Designed to maintain the size of a repaired mitral or tricuspid annulus while sustaining physiologic movement
- Anterior and posterior support allows dynamic annular motion during the cardiac cycle
- Flexible design aids cardiac output and left ventricular function, supporting optimal hemodynamics¹⁻³
- Large target area eases suturing
- Suitable for implant using open sternotomy, minimally invasive and robotic approaches

Ordering Information

Contents: Flexible Ring (1 unit per box)

Model/Reorder Number	Ring Size	Intertrigonal Dimension (mm) [A]	Inside Dimension (mm) [B]	Outside Dimension (mm) [C]
TARP-25	25	25	28.4	34.4
TARP-27	27	27	31.2	37.2
TARP-29	29	29	34.0	40.0
TARP-31	31	31	36.6	42.6
TARP-33	33	33	39.4	45.4
TARP-35	35	35	40.9	47.9



1. Yamaura Y, Yoshikawa J, Yoshida K, et al. Three-dimensional analysis of configuration and dynamics in patients with an annuloplasty ring by multiplane transesophageal echocardiography: Comparison between flexible and rigid annuloplasty rings. *J Heart Valve Dis.* 1995;4(6):618-22.
2. Kunzelman K, Reimink S, Cochran P. Flexible versus rigid ring annuloplasty for mitral valve annular dilatation: A finite element model. *J Heart Valve Dis.* 1998;7(1):108.
3. David T, Armstrong S, Sun Z. Left ventricular function after mitral valve surgery. *J Heart Valve Dis.* 1995;4(Suppl) S175-B0.

Please review the Instructions for Use prior to using these devices for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

Product referenced is approved for CE Mark.

Tailor, ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL, LESS RISK are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies.

Structural Heart International Catalogue
GMCS0177EN (May 2012)

©2012 St. Jude Medical, Inc. All Rights Reserved.

Page VR3