

1 Reducing the Risk of AV Block

Elongatable waist designed for a variety of VSD

Soft woven mesh provides high conformability to septal defects



2 Low Profile

Partnered with 4-7 Fr. delivery sheath
Slim cable to minimize unwanted damage

3 Procedural Flexibility

Double-sided screw for retrograde or antegrade approach

Ordering Information

Catalogue Number	D Disc Diameter (mm)	D2 Waist Diameter LV Side (mm)	D1 Waist Diameter RV Side (mm)	L Waist Length (mm)	Recommended Delivery Sheath
	LT-MFO-5-3	10	5	3	4F - 5F
	LT-MFO-6-4	10	6	4	
	LT-MFO-7-5	12	7	5	
	LT-MFO-8-6	12	8	6	
 with membrane	LT-MFO-9-7	14	9	7	6F
	LT-MFO-10-8	14	10	8	
	LT-MFO-12-10	16	12	10	7F
	LT-MFO-14-12	18	14	12	



LifeTech Scientific Corporation

Add: Cybio Electronic Building, Langshan 2nd Street, North Area of High-tech Park,
 Nanshan District, Shenzhen 518057, P. R. China
 Tel: +86-755-86026250 Fax: +86-755-86026251
 Email: marketing@lifetechmed.com
 Web: www.lifetechmed.com

EU authorized representative:
 LifeTech Scientific (Europe) Coöperatief U.A.
 Oudemolenstraat 60, 6416CB Heerlen, the Netherlands